


## THE **f** FRANKENMAN CK KIT PURSESTRING ANOSCOPE KIT



### CIRCULAR ANAL DILATOR 'CAD'

PRODUCT CODE

SUTURE  
THREADER

UNIVERSAL  
(FEMALE)

BUTTERFLY  
(MALE)

ANAL  
DILATOR

PURSE-STRING  
ANOSCOPE

					
CK32M	1	-	1	1	1
CK34M	1	-	1	1	1







WHEN YOU PUT PATIENTS' HEALTH FIRST,  
**REMARKABLE THINGS HAPPEN**


Frankenman were proud to be the first company to manufacture and supply a dedicated stand alone, "CK Kit" to assist surgeons in the application of a pursestring suture in complex trans-anal and low rectal procedures.

## CK KIT

Each CK Kit is available in both 32mm and 34mm diameters, to allow for variations in patient anatomy

Each CK Kit contains the following components.

-  1 x Suture Threader
-  1 x Circular Anal Dilator "CAD" – "Universal" or "Butterfly"
-  1 x Anal Dilator
-  1 x Transparent Purse-String Anoscope.




## CIRCULAR ANAL DILATOR

Frankenman were the first company to manufacture a "Circular Anal Dilator", designed specifically for patients with a narrow pelvis, where the ischium may restrict insertion, such as observed in male anatomy. Two different styles of "CAD" are available.





 Butterfly (Male)





 Universal (Female)

## PURSE-STRING ANOSCOPE


 Frankenman was the first company to recognise and realise the benefits of a **Transparent** Anoscope, in terms of improved visualisation.

 The addition of a “measuring scale” also assists in the correct placement of a pursestring suture








## USE IN TATME

 The Frankenman Purse String Anoscope Kit, is the perfect device to utilise in the application of an accurate Distal Purse string suture during the TaTME procedure.



## ASSOCIATED PRODUCTS FOR TATME

Frankenman are the only manufacturer to be able to provide custom kits, which include a full range of products for TaTME procedures, including:

-  Transanal Access Port
-  Dedicated Short, extended anvil, compact stapler
-  Anvil holder
-  Linear Cutting stapler
-  Purse-string Anoscope Kit



**REY MEDICAL SA: YOUR PARTNER IN SURGERY**

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: [www.reymedical.ch](http://www.reymedical.ch) or call us: 022 743 24 10