


SACROMESH SOFT PROLAPS®
THE SHAPE THAT SIMPLIFIES YOUR SURGERY





Sacrocolpopexy designed implant
Separate tension adjustment
Optimised patient comfort



SACROMESH SOFT PROLAPS®

Treatment of
prolapse by
sacrocolpopexy



Anterior mesh is identified
by a green thread

INDICATIONS

Sacromesh Soft Prolaps® prosthesis are designed for use as reinforcement and support of tissues when treating genital prolapse surgically by sacrocolpopexy.

They can be used in either conventional or laparoscopic surgery.

PRODUCT DESCRIPTION

Sacrocolpopexy designed implant

- Precut mesh
- Anterior mesh identification with green line

Separate tension adjustment


- 2 separate meshes for sacrocolpopexy

Optimised patient comfort

- Minimal invasive surgery
- Lightweight polypropylene
- Controlled and oriented elasticity

MATERIAL

- Knitted polypropylene monofilament





A mesh tailored for
women's anatomy

SPECIFICATIONS

Specifications	Standards
Macroporous : Size > 0.075 mm	NF S 94-801
Microporous : Size > 10 microns	NF S 94-801
Average surface mass : 39 grams /m ²	NF EN 12127
Resistance > 16 N/cm	NF EN ISO 13934-1
Average thickness : 0,38 mm	NF EN ISO 5084

REFERENCES AND SIZES

Shape	Reference	Description and sizes (cm)	Packaging
	SACROMESH9	Posterior mesh 6,4 x 19 Anterior mesh 5,4 x 19,3	1 kit/box
	GCBPRANT1U	Anterior mesh 5,4 x 19,3	1/box

Sacromesh is a class IIb medical device manufactured by Cousin Biotech S.A.S. The CE conformity has been carried out by the notified body SGS Belgium NV (CE1639). The management system of COUSIN BIOTECH S.A.S is certified for compliance with ISO 13485 standard. Please read carefully the instructions for use before using the device. The IFU is available electronically at www.cousin-biotech.com/en/implant-notice

Reference: FPGSACGB07 – 08/02/22

Non contractual pictures and texts. Specifications likely to be modified without notice.

Cousin Biotech S.A.S capital : 340 656 € - 398 460 261 RCS Lille - N°TVA FR 34 398 460 261

Cousin Surgery is a new marketing name of the company however the legal name remains Cousin Biotech.

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time