

**LAGIS®**



Safe and Effective Innovation By  
Collaborating With Surgeons

**TROCAR SYSTEM**


TROCAR SYSTEM <sup>TM</sup>  
***FineSEAL***

**REY MEDICAL SA: YOUR PARTNER IN SURGERY**

• QUALITY • SECURITY • ECONOMY

Visit our website: [www.reymedical.ch](http://www.reymedical.ch) or call us: 022 743 24 10

# LAGIS CLEARVIEW™ OPTICAL TROCAR



## ClearView™ Optical Trocar

Set include 1 Sleeve & 1 ClearView Optical Obturator

CODE	DESCRIPTION	LENGTH
TR-06V1	Stability sleeve, 5mm diameter	78mm
TR-06VL1	Stability sleeve, 5mm diameter	100mm
TRO-12V1	Stability sleeve, 12mm diameter	108mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.

# LAGIS DILATING TIP TROCAR



## Bladeless Dilating Tip

- Minimize abdominal wall and vessel trauma

*Improving Patient Outcomes*



## Dilating Tip Trocar

**Set include 1 Sleeve & 1 Dilating Obturator**

CODE	DESCRIPTION	LENGTH
TR-06D1	Stability sleeve, 5mm diameter	78mm


## Dilating Tip Trocar Kit

**Set include 2 Sleeves & 1 Dilating Obturator**

CODE	DESCRIPTION	LENGTH
TR-06D2	Stability sleeve, 5mm diameter	78mm
TR-06DL2	Stability sleeve, 5mm diameter	100mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.

# LAGIS BLUNT TIP TROCAR



## Blunt Tip Hasson Trocar

CODE	DESCRIPTION	LENGTH
TRO-12B1-AH	Smooth sleeve + Hasson Cone, 12mm	108mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.

# LAGIS UNIVERSAL TROCAR SLEEVE

## Durable Universal Seal

- No reducer caps needed




## Recessed Stopcock Valve

- Reduce unintentional opening/closing during surgery

## Funnel Lead-In

- Enable smooth instrument exchange



## True-Seal™ General Converter

- No reducer caps needed

## Flexible, Durable Seal

- Facilitate specimen removal and ensure constant insufflation

## Color-Coded Trocar Size

- No confusion between 11's & 12's

## Stability Thread Design

- Great abdominal wall retention



## UNIVERSAL TROCAR SLEEVE

CODE	DESCRIPTION	LENGTH
TR-06	Stability sleeve, 5mm diameter	78mm
TR-06L	Stability sleeve, 5mm diameter	100mm
TRO-12	Stability sleeve, 12mm diameter	108mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.