

FRANKENMAN


LINEAR STAPLER PORTFOLIO




AN UNPARALLELED BALANCE OF QUALITY AND COST EFFECTIVENESS


WHEN YOU PUT PATIENTS' HEALTH FIRST,
REMARKABLE THINGS HAPPEN

SINGLE USE RELOADABLE LINEAR STAPLER AND RELOADS




-  LARGE PARALLEL JAW OPENING - 20MM
-  CONSISTENT COMPRESSION
-  CONSISTENT STAPLE FORMATION
-  MANUAL/AUTOMATIC TISSUE LOCATION PIN
-  CARTRIDGE SAFETY INTERLOCK DEVICE
-  PREVENTS FIRING OF USED/EMPTY CARTRIDGES
-  MIDDLE INTERLOCK POSITION





- 1 RELOAD CARTRIDGE
- 2 RELEASE BUTTON
- 3 FIRING HANDLE
- 4 CLOSING HANDLE
- 5 JAWS
- 6 RETAINING PIN



ORDERING CODES AND SPECIFICATIONS

PRODUCT CODE		CLOSED STAPLE HEIGHT (mm)	OPEN STAPLE HEIGHT (mm)	STAPLE QUANTITY	STAPLE CROWN	NO. OF STAPLE ROWS	PER BOX
LS 30W*		1.0	2.5	23	3	3	3
LS 30B*		1.5	3.8	13			
LS 30G*		2.0	4.5	13			
LS 45B*		1.5	3.8	19			
LS 45G*		2.0	4.5	19			
LS 60B*		1.5	3.8	25			
LS 60G*		2.0	4.5	25			
LS 90B*		1.5	3.8	37			
LS 90G*		2.0	4.5	37			

*References available on demand

RELOAD CODE		PER BOX
LSR 30W*		12
LSR 30B*		
LSR 30G*		
LSR 45B*		
LSR 45G*		
LSR 60B*		
LSR 60G*		
LSR 90B*		
LSR 90G*		

ALSO AVAILABLE FROM FRANKENMAN

STANDARD CIRCULAR STAPLERS



COMPACT CIRCULAR STAPLERS



LINEAR STAPLERS



LINEAR CUTTING STAPLERS



TRANSCUT™ - TRANSVERSE LINEAR CUTTING STAPLERS




CIRCULAR STAPLERS FOR RECTAL PROLAPSE & ANOPEXY



UNIVERSAL TRANSANAL ACCESS PLATFORM



ENDOSCOPIC LINEAR CUTTER



ACCESSORIES




QUALITY CONTROL

We pride ourselves on world leading continual quality improvement processes.

ALL stapling devices and reloads have a unique serial number, as well a batch/lot number, to allow for more accurate tracking and traceability.

In addition we test fire 100% of our products (before they are then reloaded and packaged), to verify correct B-Shaped closure.



YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time